


## LABORATORNÍ SKLO

### BAŇKY SE ZÁBRUSEM



**Baňka na stanovení jodového čísla - se skleněnou zátkou typ 8481**

Obj. číslo	Objem [ml]	NZ [mm]	d [mm]	h [mm]	Balení [ks]
123 200 001	50	19/26	51	85	10
123 200 002	100	19/26	64	105	10
123 200 003	250	29/32	85	135	10
123 200 004	500	29/32	105	170	10
123 200 005	1000	29/32	131	215	10

**Baňka na stanovení jodového čísla - se zátkou PE**


Obj. číslo	Objem [ml]	NZ [mm]	d [mm]	h [mm]	Balení [ks]
123 200 006	50	19/26	51	85	10
123 200 007	100	19/26	64	105	10
123 200 008	250	29/32	85	135	10
123 200 009	500	29/32	105	170	10
123 200 010	1000	29/32	131	215	10



**Baňka kuželová - dle Erlenmeyera s NZ typ 8024**

Obj. číslo	Objem [ml]	NZ [mm]	d [mm]	h [mm]	Balení [ks]
------------	------------	---------	--------	--------	-------------


123 200 011	25	14/23	42	70	10
123 200 012	25	19/26	42	70	10
13 200 013	50	14/23	51	85	10
123 200 014	50	19/26	51	85	10
123 200 015	50	29/32	51	85	10
123 200 016	100	14/23	64	105	10
123 200 017	100	19/26	64	105	10
123 200 018	100	29/32	64	105	10
123 200 019	200	29/32	79	130	10
123 200 020	250	14/23	85	135	10
123 200 021	250	29/32	85	135	10
123 200 022	250	45/40	85	135	10
123 200 023	300	29/32	87	155	10
123 200 024	500	29/32	105	170	10
123 200 025	500	45/40	105	170	10
123 200 026	1000	29/32	131	215	10
123 200 027	1000	45/40	131	215	10
123 200 028	2000	29/32	166	275	6
123 200 029	2000	45/40	166	275	6



Baňka kuželová - s GL a modrým uzávěrem

Obj. číslo                      Objem [ml]                      GL [mm]                      d [mm]                      h [mm]


123 200 030	100	25	64	100
123 200 031	250	32	85	140
123 200 031	500	32	105	175
123 200 032	1000	32	131	215



Baňka s plochým dnem a NZ typ 8003

Obj. číslo	Objem [ml]	NZ [mm]	d [mm]	h [mm]	Balení [ks]
123 200 033	50	14/23	51	100	10
123 200 034	50	19/26	51	100	10
123 200 035	50	29/32	51	100	10
123 200 036	100	14/23	64	110	10
123 200 037	100	19/26	64	110	10
123 200 038	100	29/32	64	110	10
123 200 039	250	29/32	85	140	10
123 200 040	250	45/40	85	140	10
123 200 041	500	29/32	105	170	10
123 200 042	500	45/40	105	170	10
123 200 043	1000	29/32	131	200	10
123 200 044	1000	45/40	131	200	10
123 200 045	2000	29/32	166	250	6
123 200 046	2000	45/40	166	250	6
123 200 047	4000	45/40	207	300	2
123 200 048	6000	60/46	236	340	2

123 200 049	10000	60/46	279	400	1
-------------	-------	-------	-----	-----	---




Baňka s kulatým dnem a NZ typ 8004

Obj. číslo	Objem [ml]	NZ [mm]	d [mm]	h [mm]
123 200 050	50	14/23	51	105
123 200 051	50	19/26	51	105
123 200 052	50	29/32	51	105
123 200 053	100	14/23	64	115
123 200 054	100	19/26	64	115
123 200 055	100	29/32	64	115
123 200 056	250	29/32	85	145
123 200 057	250	45/40	85	145
123 200 058	500	29/32	105	175
123 200 059	500	45/40	105	175
123 200 060	1000	29/32	131	210
123 200 061	1000	45/40	131	210
123 200 062	2000	29/32	166	260
123 200 063	2000	45/40	166	260
123 200 064	4000	45/40	207	315
123 200 065	6000	60/46	236	355
123 200 066	10000	60/46	279	420




**Mikrobaňka - se zaobleným dnem a NZ typ 8026**

Obj. číslo	Objem [ml]	NZ [mm]	d [mm]	h [mm]
123 200 067	10	14/23	32	50
123 200 068	25	14/23	42	60
123 200 069	50	14/23	54	74
123 200 070	100	14/23	65	88
123 200 071	250	14/23	90	116



**Mikrobaňka srdcovitá s NZ typ 8027**


Obj. číslo	Objem [ml]	NZ [mm]	h [mm]
123 200 072	10	14/23	62
123 200 073	25	14/23	75
123 200 074	50	14/23	88
123 200 075	100	14/23	103
123 200 076	250	14/23	130



**Mikrobaňka kapkovitá s NZ typ 8029**


Obj. číslo	Objem [ml]	NZ [mm]	d [mm]	h [mm]
------------	------------	---------	--------	--------

123 200 077	10	14/23	26	60
123 200 078	25	14/23	34	80
123 200 079	50	14/23	42	100
123 200 080	100	14/23	52	123
123 200 081	250	14/23	72	170



Baňka destilační - kulaté dno, 2 postranní tubusy rovné typ 8016


Obj. číslo	Objem [ml]	NZ1 [mm]	NZ2 [mm]	NZ3 [mm]
123 200 082	250	29/32	14/23	14/23
123 200 083	500	29/32	14/23	14/23
123 200 084	1000	29/32	14/23	29/32
123 200 085	2000	29/32	14/23	29/32
123 200 086	4000	45/40	14/23	29/32
123 200 087	6000	60/46	14/23	29/32



Baňka destilační - kulaté dno, 2 postranní tubusy šikmé typ 8018


Obj. číslo	Objem [ml]	NZ1 [mm]	NZ2 [mm]	NZ3 [mm]
123 200 088	250	29/32	14/23	14/23
123 200 089	500	29/32	14/23	14/23
123 200 090	1000	29/32	14/23	29/32
123 200 091	2000	29/32	14/23	29/32

123 200 092	4000	45/40	14/23	29/32
123 200 093	6000	60/46	14/23	29/32




Baňka - s kulatým dnem, 1 postranní šikmý tubus typ 8017

Obj. číslo	Objem [ml]	NZ1 [mm]	NZ2 [mm]
123 200 094	250	29/32	14/23
123 200 095	500	29/32	14/23
123 200 096	1000	29/32	14/23
123 200 097	2000	29/32	14/23
123 200 098	4000	45/40	14/23
123 200 099	6000	60/46	14/23



Baňka sulfonační - se 3 postranními tubusy typ 8022

Obj. číslo	Objem [ml]	NZ1 [mm]	2x NZ2 [mm]	NZ3 [mm]
123 200 100	500	45/40	14/23	14/23
123 200 101	1500	71/51	29/32	14/23
123 200 102	2500	71/51	29/32	14/23
123 200 103	4000	71/51	29/32	14/23
123 200 104	6000	71/51	29/32	14/23



**Baňka kultivační - dle Blacka typ 3120**

Obj. číslo	Objem [ml]	l [mm]	l1 [mm]	h [mm]	J [mm]
123 200 105	150	45	45	155	24/20
123 200 106	200	45	45	175	24/20